NSK

Air Turbine

X450L/X450 X450KL/X450SL X450BLED/X450WLED X450WHL/X450MWL/X450STL X450M4/X450 5H

Please read this Operation Manual carefully before use and file for future reference.

MADE IN JAPAN

Symbol

This product is Autoclavable up to Max.135°C. This product can be washed via Thermo Disinfector.

C € Conforms to CE European Directive of "Medical equipment directive 93/42/EEC."

Authorised representative in the European community.

Rx Only Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

1. User and Intended Use

User: Licensed Physicians (Dentist)

Intended Use: Dental Treatment (Operations on third molars and impacted tooth removal)

2. Directions for handling and operation

- Please read these Directions carefully and use only as directed.
- Failure to follow safety instructions could result in injury to the user or patient and damage to the device. Safety instructions are classified according the the severity of risk as follows:

Classification	Severity
⚠ WARNING	A potential hazard that could result in serious injury or damage to the device if the safety instructions are not followed.
(CAUTION	A potential hazard that could result in light or moderate injury or damage to the device if the safety instructions are not followed.
NOTICE	General product information highlighted to prevent product malfunction and reduction in performance.

△ CAUTION

- · Read this Operation Manual completely before use to fully understand the product functions.
- · When operating the handpiece always consider the safety of the patient.
- · Users are responsible for the proper operation and maintenance of this product.
- Do not attempt to disassemble or tamper with the handpiece unless authorized by NSK in this Operation Manual.
- · Do not drop, strike or permit any impact to the handpiece.
- · Proper protective equipment (i.e. eyewear and masks) must be worn when operating this handpiece.
- · Should the handpiece malfunction, cease operation immediately and contact your Authorised NSK Dealer.
- · Pressing the push button while the handpiece is rotating will result in OVERHEATING of the handpiece. Contact with cheek tissue may cause the push button to be pressed and burn the patient.
- Do not use highly acidic water or sterilizing solutions to wipe, immerse or clean the handpiece.
- This product is delivered in a non-sterile condition and must be autoclaved prior to use.
- · Perform routine function and maintenance checks.
- · If the product is not used for an extended period perform a function check before use on a patient.
- · To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during use.
- · Avoid eye contact with the LED light. (X450BLED, X450WLED)

- This product is classified as a Class 1 LED Product. (X450BLED, X450WLED)
- · Use a power source which meets the following requirements. (X450BLED, X450WLED)
 - 1. The power supply should be below 15W under both normal and single-fault condition.
- 2. The power supply must use a SELV circuit.
- 3. The output voltage of the power source is within the range recommended by the manufacturer of this product.
- · U.S. Federal law restricts this device to sale by or on the order of a licensed physician.
- Under normal operating conditions air exhausts from the back of the head of the handpiece. Directing air into an open wound or beneath the mucosa or dermis may result in injury from air emphysema or embolism.
- · When using this handpiece do not direct the back of the head towards a tissue flap.

3. Setting of Air & Water Supply Pressure

Supply pressure is measured at the handpiece / hose connection point and set to the value listed on the specification table. (Fig. 1)

For Multi Gauge information refer to Option Parts List.

WARNING

Do not exceed the maximum pressure specified on the specification table.

CAUTION

Do not use air contaminated by dust, moisture or oil.

4. Connecting and Disconnecting the Handpiece

Refer to Operation Manuals of coupling and hose before connecting the handpiece.


Fig.

X450L/X450 (1) Connecting

- 1) Insert the handpiece into the Coupling.(Fig. 2)
- 2) Make sure the handpiece is fully seated in the coupling.
- (2) Disconnecting

Pull back the Retention Lock Ring and remove the handpiece from the Coupling. (Fig. 2)

Do not operate the Retention Lock Ring while pressurized. This may cause the handpiece to shoot off of the coupling unexpectedly.

X450M4/X450 5H (1) Connecting

- 1) Insert the handpiece correctly into the hose connector and tighten the hose nut. (Fig. 3)
- 2) Make sure the handpiece is securely connected to the hose.
- (2) Disconnecting
 Loosen the Hose Nut and remove from the Hose.


Fig. 2



Fig. 3

5. Insertion & Removal of the Bur

5-1 To Insert the Bur

- 1) Insert the bur until it is fully seated. (Fig. 4)
- 2) Depress the Push Button and insert the bur into the chuck until it is secure then release the button.
- 3) Ensure that the bur is secure by gently pulling and pushing the bur WITHOUT depressing the Push Button.

5-2 To Remove the Bur

Depress the Push Button firmly and remove the bur.

NOTICE

Gripping the handpiece while placing your thumb tip on the Push button will make it easier to press the button.

- · Ensure the bur is fully seated in the chuck before use.
- · Do not attempt to remove the bur until it has come to a complete stop.
 - · Do not use dirty, worn or used single use burs.
- Debris from the bur shank that is introduced into the chuck could cause the bur to slip and prevent the handpiece from fully gripping the bur.
- · Do not exceed the bur speed recommended by the bur manufacturer.
- Do not exceed maximum bur length recommended by NSK.
- Do not apply excessive lateral pressure to the bur as it may break, bend or become difficult to remove.
- Damaged burs may break at high rotation speeds or slip from the chuck. Always inspect the bur before use and do not use any bur with damage that includes but not limited to:
- · Bent, deformed, worn, rusted, broken bur.
- Cracked.
- · Non-ISO standard or tampered bur

6. Pre Check

Check that the Head Cap is securely tightened. Check for handpiece vibration, noise and overheating. If anything appears to be malfunctioning do not use the handpiece and contact your authorised NSK Dealer.

7. Maintenance

Perform the following maintenance after every patient.

7-1 Cleaning handpieces with the NSK Clean-Head System

After the treatment of each patient, clean the Clean Head.

- 1) Remove any debris from the Clean Head Holes with the provided cleaning wire and brush. (Fig. 5)
- 2) Fill a cup with clean water.
- 3) While rotating the handpiece immerse the front half of the handpiece head in the water. (Fig. 6)
- 4) Rotate the handpiece 3 times for 2 to 3 seconds each time.
- 5) Wipe the handpiece dry.

Fig. 5


Fig. 6

7-2 Cleaning (Handpiece)

- 1) Remove debris from the handpiece. Do not use a wire brush.
- 2) Wipe clean with an alcohol saturated cotton swab or cloth.
- This icon indicates that the handpiece can be washed in a Thermo Disinfector.

Refer to the Thermo-Disinfector manual for proper automatic washing procedures.

- After washing with Thermo-Disinfector and prior to lubrication ensure that the handpiece is free of any internal moisture. Moisture inside the handpiece could reduce the effect of lubrication and cause corrosion.
- · Never use and solvent such as benzine or thinner to clean the handpiece.

7-3 Cleaning (Optic)

Clean the Glass Rod tip with an alcohol-saturated cotton swab. Remove all debris and oil. (Fig. 7)

(1) CAUTION

Do not use a sharp tool to clean the Glass Rod. It could damage the glass and affect light transmission.

7-4 Cleaning (Spray Port)

If the spray ports become cloqued or cooling spray is uneven they can be cleaned using the following procedure: (Fig. 8)

\triangle CAUTION

- Do not forcibly insert the wire into the port. Resulting damage could cause the spray to be directed away from the bur and reduce the cooling effect of the spray.
- · Do not blow air into the Clean Head Hole.



Fig. 7

Fia. 8

7-5 Cleaning (Rear of the handpiece) (X450STL)

Wipe all internal surfaces at rear of the handpiece with an alcohol-saturated cotton swab.


7-6 Lubrication

NSK PANA SPRAY Plus

X450L/X450/X450KL/X450SL/X450BLED/X450WLED/X450WHL/X450MWL/X450STL

Apply NSK PANA SPRAY Plus every time after each use and/or before autoclaving.

- 1) Remove the bur from the handpiece.
- 2) Insert the Spray Nozzle onto the Pana Spray can.
- 3) Insert the Spray Nozzle in the rear of the handpiece. Hold the handpiece and spray for approximately 2-3 seconds. Continue until it runs clear from the handpiece head for at least 2 seconds. (Fig. 10)

X450M4/X450 5H

- 1) Remove the handpiece from the hose.
- 2) Insert the arrow-head spray nozzle tip into the spray can. (Fig. 11)
- 3) Place the Spray Nozzle into the Drive Air Port of the handpiece and spray for approximately 2-3 seconds. Continue spraying until it runs clear from the handpiece head for at least 2 seconds. (Fig. 12)


Fig. 11

CAUTION

- · When applying spray be sure to hold the handpiece firmly to prevent the handpiece from slipping out of the hand due to the spray pressure.
- · Hold the spray can upright.



Fig. 12

- Lubrication Oil (X450M4/X450 5H)
- 1) Remove the handpiece from the hose.
- 2) Apply 1-2 drops of lubrication oil into the Drive Air Port shown in the Fig. 12.
- 3) Connect the handpiece to the hose. Operate the handpiece for about 10 seconds.

Chuck cleaning

Clean Push Button chuck once a week.

- 1) Insert the arrow-head spray nozzle tip into the spray can nozzle.
- 2) Lubricate the chuck with NSK PANA SPRAY Plus for 1 to 2 seconds directly through the bur insertion hole. (Fig. 13)
- 3) Lubricate handpiece by using NSK PANA SPRAY Plus or Care3 Plus.

$oldsymbol{igwedge}$ Caution

If the chuck is not regularly cleaned the chuck grip may weaken and the bur may slip during use.

NSK Care3 Plus

Refer to Care 3 Plus system instructions for proper use.

7-7 Sterilization

The only recommended sterilization method is by steam autoclaving. Remove the bur after each patient and sterilize as follows.

- 1) Place the handpiece into an autoclave pouch and seal.
- Autoclave according to the following time requirements.
 Autoclave for more than 20 min. at 121°C. 15 min. at 132°C or 3 min. at 134°C.
- 3) The handpiece should remain sealed in the autoclave pouch until required for use.

CAUTION

- Do not autoclave the handpiece with other instruments even in a pouch. This is to prevent possible discoloration and damage to the handpiece from chemical residue on other instruments
- · Keep the handpiece in normal atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- · Immediately after use the handpiece should be cleaned, lubricated and sterilized. If blood remains on any external or internal surfaces it can harden and cause rust.
- · Do not heat or cool the handpiece too quickly. Rapid change in temperature could cause damage to the handpiece.

- Do not use the dry cycle if the sterilizer chamber temperature may exceed 135°C.
- · Autoclave sterilization is recommended for this product. The validity of other sterilization methods is unknown.
- · Do not touch the handpiece immediately after autoclaving as it will be very hot.


NOTICE

NSK recommends Class B sterilizers as stated in EN13060.

8. Replacing the 0-rings (X450L/X450)

Replace the O-rings if water is present in the exhaust air line, ALWAYS change the complete set of O-rings.

1) Loosen and Remove the taper ring at the rear of the handpiece. (Fig. 14)


- 2) Gently remove each O-ring by hand. (Fig. 15)
- 3) Insert the complete set of new O-rings in the correct grooves.
- 4) Replace and firmly tighten the taper ring.
- *Refer to Spare Parts List to identify the correct parts.

CAUTION

- · Do not force the new replacement O-ring with excessive pressure.
- · When inserting new O-rings, make sure they are inserted in the correct grooves.
- · Make certain that the taper ring is firmly tightened. If the taper ring is loose water and air leakage could occur.

9. Replacing the Non-Retraction Valve

A water Non-Retraction Valve is integrated in the Coupling Joint (the product below), which shuts off the water retraction directly at the handpiece head to prevent fluids infiltrating the water line. If water is beginning to leak from the handpiece, replace the non-retraction valve.

NSK Coupling

- 1) Remove the Coupling Joint from the hose.
- 2) Remove the back end Gasket.
- 3) Pull and remove the water tube, and replace the Non-Retraction Valve.
- 4) Insert the new Non-Retraction Valve securely and remount the back-end Gasket.
- *Refer to Spare Parts List to identify the correct parts.

X450BLED

- 1) Remove the handpiece from the coupling.
- 2) Loosen and remove the Taper Ring.
- 3) Pull out the joint.
- Pull out the old Non-Retraction Valve with a thin instrument. Insert a new Non-Retraction Valve securely.
- 5) Insert the joint by aligning the joint holes with the pipes of the handpiece.
- 6) Firmly tighten the Taper Ring while lightly pressing the bulb holder by a finger.
- *Refer to Spare Parts List to identify the correct parts.


Fig. 16


Fig. 17


Fig. 18

- Tighten the taper ring firmly. A loose fitting taper ring could result in water and air leakage.
- · Do not remove the internal circuit board when removing the taper ring.

10. Routine Maintenance Checks

Perform routine maintenance checks every three months. If any abnormalities are found, contact your authorized NSK Dealer.

Points to check	Details			
Head cap is loose	Check that the Head Cap is firmly tightened.			
Rotation Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.				

11. Specifications

Model	X450L	X450	X450KL/X450SL	X450BLED/X450WLED	X450WHL/X450MWL X450STL	_ X450M4	X450 5H	
Hose Connection Type	-				IS0916	68 Type B		
Max. Rotation Speed		380,000-450,000min ⁻¹						
Bur Type		ISO 1797-1 ø1.59-1.60mm Surgical Long Shank Bur (25mm)						
Chucking Length		10.7mm						
Max. Bur Length		25mm						
Max. Working Part Diameter		ø1mm						
Drive Air Pressure							0.18-0.22MPa (1.8-2.2kgf/cm²) (26.1-31.9psi)	
Max. Air Consumption	45±5NL/min (1.41±0			0.18cfm) (0.25MPa)		n 35±5NL/min (1.2	3±0.18cfm) (0.22MPa)	
Water Pressure		0.1-0.2MPa (1.0-2.0kgf/cm ²)			n ²) (14.5-29.0psi)			
Optic	Glass Rod	-	Glass Rod	White LED	Glass Rod	-	Glass Rod	
Voltage	oltage - AC/DC 3.3 ±0.05V -		-					
Active Voltage	-		2.8 - 4.0V	-				
Current	-		0.38A Typical(3.3V)	al(3.3V) -				
Use Environment	Temperature: 10-40°C			Humidity: 30-75% A	tmospheric Pressure: 700-1,060	0hPa		
Transportation and storage Environment								
*This boundaines in a contaction to the con-	*This handsigns is a water ist type. The coolant water is not approved.							

^{*}This handpiece is a water jet type. The coolant water is not sprayed.

12. Warranty

NSK products are warranted against defects in materials and workmanship. NSK reserves the right to inspect all product to determine the cause of any malfunction. Warranty is voided should the handpiece be used incorrectly, for anything other than the intended purpose, has been tampered with by unqualified personnel or has had non-genuine NSK parts installed. Replacement parts will remain available for seven years following discontinuation of the model.

13. Option Parts List

Model	Order Code
MG-4H Multi Gauge	Z109400

14. Spare Parts List

Model	Order Code	Compatible Product
PTL O-ring Set	Y900580	X450L/X450
Non-Retraction Valve	P401054	NSK coupling
Non-Retraction Valve	P826B042	X450BLFD

15. Disposing product

In order to avoid health risks to operators handling the disposal of medical equipment as well as the risk of environmental contamination, a surgeon or a dentist is required to confirm the equipment is sterile. Seek the assistance of specialist firms who are licensed to dispose of specially controlled industrial and medical wastes to dispose the product for you.

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NSK France SAS www nsk fr
NSK Europe GmbH □ □
NSK United Kingdom Ltd
NSK Oceania Pty Ltd
NSK America Corp
NSK Asia Pte Ltd